


# JuLI™ Stage

Real-time live cell imaging system

**See the whole story,  
Do not miss a moment of your cells.**



# See what happened inside an incubator with JuLI™ Stage




**“You can do experiments  
without any limitation  
inside your incubator!”**

- Compact and compatible with a standard CO<sub>2</sub> incubator
- Fully automated X-Y-Z stage
- Multi-channel fluorescence imaging
- Easy & powerful software
- Take and analyze images in real time







## What you can do with JuLI Stage




**Spheroids**




**Apoptosis**




**Cellular localization**




**Neurite growth**




**Cytotoxicity**




**Scratch assay**




**Fluorescent cell counting**




**Proliferation**




**Transfection efficiency**




**iPS cell line**




**Embryonic morphology**



**Reporter gene**



**Chemotaxis**



**Oocyte monitoring**



**Tissue observation  
(auto-stitched images)**

# CHR

Cell History Recorder

**JuLI™ Stage, the new standard of CHR (Cell History Recorder) is designed to get time-lapse images.**


- Observe cells in real time and record a cell history from the beginning until the end
- Revert to the time point you desire
- Save time with a fully automatic time-lapse imaging function



- Well plates (6 to 384 wells)
- Slides
- Petri dishes & flasks


Compatible with various brand well plates with the auto adjustment function

## 01 Multi-Position




Take any number of images of any positions of a well

## 02 Stitching options



Easy to acquire a whole well stitching image with the stitching function

## 03 Focus options



Acquire high-quality images from the Z-stack focus option

# EDIT

Image Edit

JuLI™ EDIT is designed to import, verify and re-edit project data taken from JuLI™ Stage.

- Edit images in projects
- Make movies
- Review data on your personal PC

## 01 Image Editor




Image. Background Correction

### [ Bright field ]

Improve image quality with background correction




Image. Auto Adjust

### [ Fluorescence field ]

Improve image quality with the auto adjust function

## 02 Movie Maker



Make various types of movie


# STAT

Image Statistics

JuLI™ STAT is designed to import, verify and analyze project data taken from JuLI™ Stage.


- Get result data from projects
- Make graphs using result data

## 01 Attached Cell Counting



Fluorescent cell counting in real time

## 02 Quantitative Results



Analyze results in various forms


## 03 Plate Editor

|   | 01   | 02   | 03   | 04 | 05   | 06   | 07 | 08 | 09   | 10   | 11 | 12 |
|---|--|--|--|----|--|--|----|----|--|--|----|----|
| A | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    |
| B | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |  |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    |
| C | DMF 100.00 µM<br>Media (1) 20.00 µM / well |  |  |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    |
| D | DMF 100.00 µM<br>Media (1) 20.00 µM / well |  |  |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    |
| E | DMF 100.00 µM<br>Media (1) 20.00 µM / well |  |  |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    |
| F | DMF 100.00 µM<br>Media (1) 20.00 µM / well |  |  |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    |
| G | DMF 100.00 µM<br>Media (1) 20.00 µM / well |  |  |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    |
| H | DMF 100.00 µM<br>Media (1) 20.00 µM / well |  |  |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    | DMF 100.00 µM<br>Media (1) 20.00 µM / well | DMF 100.00 µM<br>Media (1) 20.00 µM / well |    |    |

Design your experimental workflow

## SCRATCH STAT

Scratcher easily creates a uniform scratch line for 96 wells. The Scratch STAT software analyzes scratches in real time.






[ Before : 0 hr ]

[ After : 24 hr ]

[ Scratch graph of 96 well plate ]

## SPHEROID STAT

The Spheroid STAT software provides various analysis functions of up to 96 spheroids in real time.



※ For the single spheroid cell only

## Specifications


| Items                 | Specification  |
|-----------------------|--|
| Light source          | Blue, Green, UV LED (Intensity adjustable)   |
| Objective lens        | 4X, 10X, 20X, + Digital zoom<br>Inter-changeable objective lens  |
| Fluorescence          | 3 fluorescences<br>DAPI: Excitation 378/52, Emission 447/60<br>GFP : Excitation 466/40, Emission 525/50<br>RFP : Excitation 543/22, Emission 580LP   |
| Camera                | High-sensitivity monochrome CCD (Sony sensor 2/3")<br>1,936 x 1,456 pixels (2.8 M), 53 FPS, 14 bit   |
| Stage                 | Automated, motorized, X-Y-Z stage<br>Ex-changeable vessel holders(optional)  |
| Exported formats      | Image : JPEG, TIFF, BMP, PNG<br>Video : AVI<br>Raw data : CSV  |
| PC                    | Desktop computer, Desktop monitor 24-in. LCD<br><br>CPU: Intel i5, 9 generation or over spec.<br>OS: Windows® 10 Pro 64 bit<br>RAM: 16 GB<br>Hard drive: 2 TB<br>Network: Gigabit Ethernet, WiFi<br><br><i>*PC specifications may change without notice.</i> |
| Operating power       | 100 - 240 V, 1.5 A, 50/60 Hz   |
| Electronic input      | 12 VDC, 2.0 A  |
| Operating environment | 5 - 40 °C, 20 - 95%  |
| Dimensions            | 429(W) X 310(D) X 324(H) mm  |
| Weight                | 18.5 kg / 41 lb  |

## Ordering Information

| Cat. No. | Product   | Description  |
|----------|---|--|
| JS1000S  | JuLI™ Stage, Starter Pack                       | JuLI™ Stage basic set (JS1000), Desktop computer (JP0200), 3 Objective lenses (4X, 10X & 20X)  |
| JS1000   | JuLI™ Stage, Real-Time Live Cell Imaging System | Main device, power supply, control box   |
| JP0100   | Desktop Computer                                | CPU: Intel i5, 9 generation or over spec.<br>OS: Windows® 10 Pro 64 bit<br>RAM: 16 GB<br>Hard drive: 2 TB<br>Network: Gigabit Ethernet |
| JMO100   | Desktop Monitor                                 | 24" Full HD (1920 x 1080) monitor  |
| JP0150   | External Hard Disk Drive (Optional)             | Total 8 TB (4 TB x 2 ea)   |
| JO0004   | Objective Lens (4X)                             | Magnification : 4X, NA : 0.16  |
| JO0010   | Objective Lens (10X)                            | Magnification : 10X, NA : 0.3  |
| JO0020   | Objective Lens (20X)                            | Magnification : 20 X, NA : 0.45  |
| JVH001   | Vessel Holder (Optional)                        | Micro Slide (26 x 76 mm)   |
| JVH002   | Vessel Holder (Optional)                        | Petri Dish (35 mm)   |
| JVH003   | Vessel Holder (Optional)                        | Petri Dish (60 mm)   |
| JVH004   | Vessel Holder (Optional)                        | Petri Dish (100 mm)  |
| JVH005   | Vessel Holder (Optional)                        | T-Flask (25 & 75 cm <sup>2</sup> )   |
| JSCT100  | JuLI Analysis Software (Scratch)                | JuLI Scratch STAT<br>JuLI Scratcher  |
| JSPT100  | JuLI Analysis Software (Spheroid)               | JuLI Spheroid STAT   |

# JuLI™ Stage

Real-time live cell imaging system



website | [www.nanoentek.com](http://www.nanoentek.com)  
e-mail | [sales@nanoentek.com](mailto:sales@nanoentek.com)

FOR RESEARCH USE ONLY.  
This product is not approved for diagnostic or therapeutic use.

## NanoEntek, Inc.

### Head Office

12F, 5, Digital-ro 26-gil, Guro-gu, Seoul, 08389, Korea  
Tel +82-2-6220-7940 / Fax +82-2-6220-7999

## NanoEntek America, Inc.

220 Bear Hill Road, Suite 102, Waltham, MA 02451, USA  
Tel +1-781-472-2558 / Fax +1-781-790-5649

## NanoEntek Europe | med-tech supplies GmbH

Lochamerstr. 4a, 82152 Martinsried, Germany  
Tel +49-89-21-55-38-43 / Fax +49-89-99-95-46-60

NESCT-JST-001E (V.2.1)