

Octet instruments have been widely deployed in various vaccine development research programs (HIV^{1,2}, influenza³, Ebola⁴, novel coronavirus⁵⁻¹⁵) and in evaluating vaccine agents such as nanoparticles¹⁶, virus-like particles (VLPs)¹⁷ and bispecific antibodies¹⁸. Octet systems also provide solutions beyond lead candidate selection to bioprocess development and manufacturing. This review summarizes the use of Octet assays in some recent advances in coronavirus research and vaccine development efforts.

Determination of coronavirus receptor binding mechanisms and cross-reactivity

The 2019-nCoV or SARS CoV-2 shares genetic and morphologic features with other coronavirus families, particularly from the betacoronavirus genus (SARS-CoV 2003). SARS CoV-2 makes use of a glycosylated, homotrimeric class I fusion spike (S) protein to gain entry into host cells via the human angiotensin converting enzyme-2 (ACE2) receptor¹⁹. The SARS CoV-2 S glycoprotein bears significant structural homology with SARS-CoV compared to other major coronaviruses such as the MERS-CoV. The S-protein is known to exist in a metastable pre-fusion conformation that undergoes rearrangements to fuse the viral membrane with the host cell membrane. Binding of the S1 subunit to host-cell receptors triggers destabilization of the pre-fusion trimer and conformational change in the protein initiates fusion with the host cell membrane through the S2 subunit. In order to engage the host-cell receptor, the receptor-binding domain (RBD) of S1 would need to undergo hinge-like transitions that either hides or displays the receptor-binding regions which are termed as *down* and *up* conformations respectively, providing accessibility to receptor-binding and presenting itself as a key antigenic target for vaccine development.

One of the first CryoEM structures of the SARS CoV-2 S trimer in the pre-fusion conformation was revealed by Wrapp *et al.* where they showed that while the overall structure of SARS CoV-2 S protein was mostly similar to SARS-CoV S, differences existed between the down conformations at the RBD¹².

Another cryo-EM structure of SARS CoV-2 spike glycoprotein complex resolved by Walls *et al.* showed that that the SARS-CoV-2 S glycoprotein harbors a furin cleavage site at the boundary between the S1/S2 subunits, which is processed during biogenesis and differentiates this virus apart from SARS-CoV and other SARS-related CoVs¹⁰. Using Octet binding assays, they showed that the receptor-binding domains of SARS-CoV-2 S and SARS-CoV S engaged the human ACE2 with similar binding affinities. Similar observations were reported by Joyce *et al.* using a high-resolution crystal structure of the SARS-CoV-2 S RBD⁵. Additionally, in this work antibodies that were known to interact with SARS-CoV and MERS CoV receptor binding domains (RBD) were tested for SARS-CoV-2 binding activity using an Octet kinetic assay. Only two antibodies, 240CD and CR3022 out of the pool, displayed low nanomolar binding affinities against the SARS-CoV-2 RBD. These two antibodies, however, competed for the same SARS-CoV-2 RBD binding site according to an Octet cross-com-

Epitope binning formats and assay guidance on Octet instrument




Figure 1: Assay formats of epitope binning assays and recommendation guide to select the most suitable epitope binning assay format. Octet systems enable multiple assay formats for cross-competition assessment and provide dedicated high-throughput software analysis tools for epitope binning. (A) In Tandem Assay: Antigen is immobilized on the biosensor, followed by the binding of the saturating mAb (Ab1) and competing mAb (Ab2), respectively. (B) Classical Sandwich Assay: One mAb is immobilized on the biosensor (Ab1), antigen is captured using this mAb and then a second sandwiching mAb is tested (Ab2). (C) Premix Assay: One mAb is immobilized on the biosensor (Ab1), the biosensor is then exposed to a premix solution containing the antigen and a large molar excess of the second mAb (Ab2). (D) Guidance for selecting the most suitable epitope binning assay format. Please read the [Octet assay guide](#) to learn more about developing and analyzing cross-competition screens.

petition assay. Cross-competition or epitope binning assays can be used to segment monoclonal antibodies (mAbs) into bins based upon the antigen region or epitope bound by each antibody (Figure 1). Furthermore, an Octet competition assay revealed that 240CD or CR3022 antibody binding to SARS CoV-2 RBD did not perturb ACE2 recognition. Crystallography data revealed that CR3022 engaged the RBD at a site that is conserved in both SARS CoV and SARS CoV-2 explaining the cross-reactivity. However, a comparison of SARS-CoV and MERS-CoV epitope regions from previous structural data to SARS CoV-2 identified in this work revealed that CR3022 was binding to a novel epitope within SARS CoV-2 S. A specific RBD knockout mutant intro-

duced within this epitope region (glycan sequon at position 384) eliminated binding to both CR3022 and 240CD antibodies according to Octet data confirming shared epitopes between both antibodies.

To understand the SARS CoV-2 spike conformations with respect to receptor binding and cell entry, structures of CR3022 binding to trimeric units of SARS-CoV-2, SARS-CoV, and MERS-CoV were modeled⁵. Data indicated that CR3022 epitopes were obstructed by adjacent spike promoters when the RBD is in the *down* conformation but were more accessible in the *up* confirmation. Octet assays observed CR3022 binding to non-stabilized S-glycoprotein (S1) conformations and much weaker binding to S2 P trimers. To assess whether minimal proteolytic action or receptor binding could increase the availability of the obscure CR3022 epitope, S2 trimers were treated with trypsin or incubated with ACE2. Incubation of S2 P trimer with human ACE2 did not dramatically affect CR3022 binding. However, trypsin treatment of the S2 P conformation resulted in increased binding levels similar to non-stabilized S glycoprotein binding, and the level of binding was titratable with increasing proteolytic action, inferring to the cryptic nature of the CR3022 binding epitope. Similar observations were reported by Yuan *et al.*, using the crystal structure of neutralizing antibody CR3022 bound to SARS-CoV-2 RBD¹⁴. Octet binding data showed that despite high sequence conservation between SARS-CoV and SARS CoV 2 RBDs, CR3022 FABs binds SARS-CoV with ~100-fold higher affinity than to SARS CoV-2 RBD possibly driven by the non-conserved residues between the two epitopes mentioned above. An *in vitro* microneutralization assay indicated higher neutralization activity with SARS-CoV, but not SARS-CoV-2 for CR3022, a result that is consistent with lower RBD binding affinities. Overall, these structural studies provide insight into how SARS CoV-2 can be targeted by the humoral immune response and revealed a conserved, but cryptic epitope shared between SARS CoV-2 and SARS CoV that can potentially be utilized to create cross-protective nAbs.

Neutralization antibody development

Neutralizing antibodies (nAbs) are likely to be one of the most effective treatments against coronavirus infections among the current therapeutic options and are in urgent need. Several nAbs targeting SARS-CoV RBD, have been shown to exhibit significant *in vivo* antiviral activity by reducing virus titers in animal models. However, the availability of nAbs with cross-reactivity could be key for sustained utilization as therapeutics that can stand against virus mutations that are typical of single standard RNA viruses.

Due to the ~80% sequence similarity between the two Spike (S) proteins of SARS CoV and SARS CoV-2, understanding whether antibodies raised from SARS-CoV spike (S) protein immunization would retain cross-reactivity to the new SARS-CoV-2 will offer important insights and guidance to therapeutic antibody and prophylactic vaccine development.


Tian *et al.*, in a pioneering study, expressed and purified the CoV-2 RBD and tested the binding of several SARS-CoV specific neutralizing antibodies that were known to target the RBD and exert potent neutralization activities against SARS-CoV⁹. SARS-CoV nAbs, m396, CR3014, CR3022 and the MERS-CoV nAb m336 were tested for SARS-CoV-2 RBD recognition (Figure 2). While most antibodies did not show appreciable binding to the SARS CoV-2 RBD, CR3022 that was isolated from a SARS CoV patient bound with a 6.3 nM affinity (k_{on} of $1.84 \times 10^5 \text{ Ms}^{-1}$ and k_{off} of $1.16 \times 10^{-3} \text{ s}^{-1}$).

Octet BLI competition assays showed that CR3022 did not compete for the ACE2 binding site in SARS-CoV-2 RBD, indicating binding to an epitope different from the ACE2 binding site (Figure 2B). Similarly, Wrapp *et al.* investigated another set of previously published SARS-CoV RBD binding antibodies, S230, m396, and 80R with SARS CoV 2 S protein¹². Octet binding assays indicated a lack of significant binding activity possibly due to the recognition through different epitopes compared to SARS CoV. Sun *et al.* reported that several pAbs and mAbs generated by immunizing mouse and rabbits with SARS-CoV S1 or RBD proteins did not show favorable cross-neutralization activities in SARS-CoV-2 pseudoviruses (PSV). However, these antibodies exhibited potent neutralization potencies against SARS CoV virus and high-affinity binding to SARS CoV RBD ($K_D \sim 9\text{--}200 \text{ pM}$)⁸. These results indicate that antibodies that exhibit SARS-CoV neutralization potencies or RBD binding is not necessarily a prerequisite for activity in SARS CoV-2 despite high sequence homologies at the RBD and suggests the necessity to develop novel antibodies that specifically target SARS CoV-2 RBDs.

Wang *et al.* reported the discovery of a human monoclonal antibody with neutralization potencies for SARS-CoV-2 and SARS-CoV virus activity¹¹. 47D11, derived from a screen that consisted of a collection of supernatants that contained SARS S hybridoma derived from immunized transgenic mice, showed cross-reactivity and neutralization activity against both SARS CoV and SARS CoV-2 pseudotyped VSV infections. Octet studies showed that 47D11 binds the SARS CoV S_{ecto} domain with higher affinity ($K_D = 0.745 \text{ nM}$) relative to SARS CoV-2 S_{ecto} domain ($K_D = 10.8 \text{ nM}$) although similar binding affinities were reported for binding to the SARS CoV and SARS CoV-2 RBD binding domains, further validating the differences in epitope accessibility between the two virus receptor binding motif (RBM) conformations. Additionally, as in other cases reported thus far, 47D11 did not compete with the SARS CoV or SARS CoV-2 RBD, implying neutralization through a mode that is different from receptor blocking. This cross-neutralizing antibody reportedly targets a communal epitope and offers potential for prevention and treatment of SARS CoV-2.

Pinto *et al.* screened for nAbs from a SARS CoV survivor⁷. Eight mAbs from the screen bound to both SARS CoV and SARS CoV-2 S protein-transfected CHO cells, out of which mAbs S303, S304, S309, and S315 recognized both SARS CoV and SARS CoV-2 RBDs. Octet data indicated S309 binding to both S domains with similar nanomolar affinities and also showed

A



B




Figure 2: Characterizing of SARS CoV antibody interactions to SARS CoV-2 RBD. (A) Binding profiles of SARS CoV-2 RBD to ACE2 and antibodies and (B) competition of CR3022 and ACE2 with SARS CoV-2 RBD measured by Octet. For the competition analysis, ACE2 were immobilized onto HIS1K biosensors, followed by binding to either ACE2, CR3022 + ACE2 or isotype antibody control of ACE2. Figure is reproduced from reference 9.

comparable neutralization potencies against both SARS CoV and CoV-2 pseudoviruses and authentic SARS-CoV-2. According to structural data, S309 recognizes a protein/glycan epitope on the SARS CoV-2 RBD distinct from the receptor-binding motif (RBM) and seems to be accessible to both up and down states. Octet BLI competition assays showed that binding of S309 Fab and IgG forms to SARS CoV-2 RBD did not discourage ACE2 receptor binding to SARS CoV-2 S protein. To gain more insight on the epitopes recognized by the panel of mAbs, an Octet epitope binning assay was carried out to map the antigenic sites present on the SARS CoV and, SARS-CoV-2 RBDs. Octet epitope binning analysis identified four antigenic regions within the RBD of SARS CoV where the mAbs were targeting. Based on these findings, mAb combinations that targeted different antigenic regions were tested for neutralization potencies and possible synergistic effects where data indicated that either S304 or S315 in combination with S309 enhanced neutralization potencies against SARS CoV-2. S309 identified in this work showed broad neutralization activity against multiple sarbecoviruses showing promise as an effective therapeutic candidate.

Pan *et al.* used the SARS CoV-2 RBD domain to produce neutralization antibodies in horse antisera⁶. F(ab')₂s isolated from horse antisera reportedly showed neutralization activity against SARS CoV-2 with an EC₈₀ of ~25 µg/mL. Octet BLI kinetic analy-

sis also confirmed SARS CoV-2 receptor binding to F(ab')₂ with an affinity of 76 nM. Affinity purification of F(ab')₂ from antisera improved neutralization activity against SARS CoV-2 (EC₈₀ = 0.18 µg/mL) and increased binding affinity (0.76 nM) measured by Octet BLI, highlighting RBD-specific F(ab')₂ as a potential therapeutic candidates for SARS CoV-2.

Wu *et al.* reported the generation and testing of single-domain antibodies targeting SARS CoV-2 RBD¹³. The panning using SARS-CoV-2 RBD and S1 as antigens resulted in the identification of antibodies targeting five types of neutralizing or non-neutralizing epitopes on SARS-CoV-2 RBD. They tested 18 human single-domain antibodies in competition binding assays using the Octet system and found that they could be divided into three competition groups (groups A, B, or C) that did not show any competition with each other. The single domain antibody n3130 identified from the library using SARS-CoV-2 S1 as the panning antigen showed potent neutralization activity with both pseudotyped and live virus and strong binding to SARS CoV-2 S1 and RBD domains. This work also reports the unique immunogenic profiles of SARS-CoV-2 RBD compared to that of SARS-CoV and MERS-CoV, which may have important implications for the development of effective vaccines against SARS CoV-2.

Design and validation of a peptide inhibitor to block viral entry

In addition to protective antibodies, peptide and small molecule inhibitors can also be used as viable inhibitory strategies to target the SARS CoV-2 RBD to prevent virus entry. Zhang *et al.* investigated the design and development of a 23-mer peptide inhibitor (SBP1) discovered based on the ACE2 PD α1 helix that spans the binding interface according to the crystallographic structure of the SARS-CoV-2 S receptor–ACE2 protein complex¹⁵. Octet BLI was utilized to test target specificity of the peptide binding to SARS CoV-2 RBD. SBP1 bound the RBD with an affinity of 47 nM ($k_{on} = 4.69 \times 10^4 \text{ M}^{-1} \text{ s}^{-1}$ and $k_{off} = 2.2 \times 10^{-3} \text{ s}^{-1}$) where this binding affinity is comparable to that of the full-length ACE2 binding to SARS CoV-2 RBD. SBP1 could, therefore, potentially cover spike proteins on the SARS-CoV-2 surface and outcompete the binding for ACE2 as a novel therapeutic intervention strategy.

Summary

Affinity characterization, neutralization, and competition studies are essential in the selection of high affinity and potentially neutralizing lead candidates for therapeutic and vaccine development against the novel SARS CoV-2. Kinetic analysis further describes the components of association and dissociation that comprise the overall affinity, while cross-competition assays allow for the identification of more optimal candidates. In addition to the studies illustrated above, researchers have also utilized Octet assays to characterize and investigate other members of the coronavirus or associated species where the observations and results are invaluable towards implementing and framing current therapeutic strategies for SARS CoV-2. Table 2 provides a summary of selected studies illustrating the use of Octet assays on current and historical coronavirus research²⁰⁻²⁷. In all these studies, the Octet system's ease-of-use combined with the availability of a variety of ready-to-use biosensor chemistries has enabled users to significantly reduced time to results and advance the continued search for vaccines and therapeutics for the current pandemic.

Analyte(s)	Immobilized ligand(s)	Virus species	Assay focus	Biosensor	Reference
ACE2	Fc-tagged SARS CoV-2 RBD	SARS CoV-2	Cross-reactivity and receptor binding	AHC	12
ACE2	SARS-CoV-2 S _B SARS-CoV S _B	SARS CoV-2	Receptor binding	HIS1K	10
mAb, MERS S protein, VHH	SARS CoV-2 S1, SARS CoV-2 RBD	SARS CoV-2	VHH and mAb binding assessment, Epitope binning	AR2G	13
Fab, IgG	SARS CoV-2 S protein, SARS CoV-2 RBS	SARS CoV-2	FAB, IgG binding characterization	HIS1K	14
ACE2, SARS-CoV RBD	SARS-CoV-specific neutralizing antibodies	SARS CoV-2	Cross-reactivity binding	SA	9
Antibodies	Biotin-Recombinant SARS-CoV S1 protein	SARS CoV-2	Neutralization antibody assessment	SA	8
FAB, ACE2, SARS S protein	mAb, ACE2	SARS CoV-2	Receptor binding, antibody reactivity, competition assays	AHC, HIS1K	5
SARS CoV-2 RBD	Inhibitor Peptide	SARS CoV-2	Peptide inhibitor development	SA	15
mAb	S1B and SARS S _{ecto} domains	SARS CoV and SARS CoV-2	Antibody reactivity assessment	HIS1K	11
SARS CoV-2 RBD, SARS CoV RBD, mAbs, FAB	mAbs, SARS CoV-2 RBD, SARS CoV RBD	SARS CoV and SARS CoV-2	Antibody binding characterization, epitope binning, competition binding	Protein A, HIS2, HIS1K	7
Coronavirus S proteins	Biotin-9OAc6SLN	SARS CoV and MERS CoV	Receptor binding	SA	20
MERS-5HB fusion inhibitor	Biotin-peptide (MERS-HR2P)	MERS CoV	Protein Inhibitor development	SA	21
mAb, MERS S protein	mAb, MERS S protein	MERS CoV	Epitope binning and Receptor binding	SAX, Protein A	22
MERS-CoV S	mAbs	MERS CoV	Antibody reactivity measurements, competition assays	AHC, HIS1K	23
MERS-CoV NTD, mutants	mAb	MERS CoV	Antibody reactivity measurements	AHC	24
MERS S-protein	mAbs	MERS CoV	Antibody reactivity measurements	AHC	25
RBD/S1/S2 antigen	Vaccine-induced mouse monoclonal IgGs	MERS CoV	Antibody reactivity and epitope binning	AHC, HIS1K	26
FABs	MERS-CoV RBD	MERS CoV	Epitope binning	HIS1K	27

Table 2: Examples of Octet BLI assays in coronavirus research and vaccine development.

References

- 1 Structures of protective antibodies reveal sites of vulnerability on Ebola virus, Murin CD et al., *Proc Natl Acad Sci USA*, 111:17182-17187, 2014, doi:10.1073/pnas.1414164111.
- 2 A next-generation cleaved, soluble HIV-1 Env trimer, BG505 SOSIP.664 gp140, expresses multiple epitopes for broadly neutralizing but not non-neutralizing antibodies, Sanders RW, et al., *PLoS Pathog*, 9:e1003618, 2013, doi:10.1371/journal.ppat.1003618.
- 3 Kinetic analysis of the influenza A virus HA/NA balance reveals contribution of NA to virus-receptor binding and NA-dependent rolling on receptor-containing surfaces, Guo H, et al., *PLoS Pathog*, 14:e1007233, 2018, doi:10.1371/journal.ppat.1007233.
- 4 Antibody Repertoires to the Same Ebola Vaccine Antigen Are Differentially Affected by Vaccine Vectors, Meyer M, et al., *Cell Rep*, 24:1816-1829, 2018, doi:10.1016/j.celrep.2018.07.044.
- 5 A Cryptic Site of Vulnerability on the Receptor Binding Domain of the SARS-CoV-2 Spike Glycoprotein, Joyce MG, et al., *bioRxiv*, 2020.2003.2015.992883, 2020, doi:10.1101/2020.03.15.992883.
- 6 Immunoglobulin fragment F(ab)₂ against RBD potently neutralizes SARS-CoV-2 *in vitro*, Pan X, et al., *bioRxiv*, 2020.2004.2007.029884, 2020, doi:10.1101/2020.04.07.029884.
- 7 Structural and functional analysis of a potent sarbecovirus neutralizing antibody, Pinto D, et al., *bioRxiv*, 2020.2004.2007.023903, 2020, doi:10.1101/2020.04.07.023903.
- 8 SARS-CoV-2 and SARS-CoV spike-RBD structure and receptor binding comparison and potential implications on neutralizing antibody and vaccine development, Sun C, et al., *bioRxiv*, 2020.2002.2016.951723, 2020, doi:10.1101/2020.02.16.951723.
- 9 Potent binding of 2019 novel coronavirus spike protein by a SARS coronavirus-specific human monoclonal antibody, Tian X, et al., *Emerg Microbes Infect*, 9:382-385, 2020, doi:10.1080/22221751.2020.1729069.
- 10 Structure, function, and antigenicity of the SARS-CoV-2 spike glycoprotein, Walls, AC et al., *Cell*, 181:281-292 e286, 2020, doi:10.1016/j.cell.2020.02.058.
- 11 A human monoclonal antibody blocking SARS-CoV-2 infection, Wang C, et al., *bioRxiv*, 2020.2003.2011.987958, 2020, doi:10.1101/2020.03.11.987958.
- 12 Cryo-EM structure of the 2019-nCoV spike in the prefusion conformation, Wrapp D, et al., *Science*, 367:1260-1263, 2020, doi:10.1126/science.abb2507.
- 13 Fully human single-domain antibodies against SARS-CoV-2, Wu Y, et al., *bioRxiv*, 2020.2003.2030.015990, 2020, doi:10.1101/2020.03.30.015990.
- 14 A highly conserved cryptic epitope in the receptor-binding domains of SARS-CoV-2 and SARS-CoV, *Science*, doi:10.1126/science.abb7269 (2020).
- 15 The first-in-class peptide binder to the SARS-CoV-2 spike protein, Zhang G, et al., *bioRxiv*, 2020.2003.2019.999318, 2020, doi:10.1101/2020.03.19.999318.
- 16 Presenting native-like trimeric HIV-1 antigens with self-assembling nanoparticles, He L, et al., *Nat Commun*, 7:12041, 2016, doi:10.1038/ncomms12041.
- 17 Exposure of epitope residues on the outer face of the chikungunya virus envelope trimer determines antibody neutralizing efficacy, Fong RH, et al., *J Virol*, 88:14364-14379, 2014, doi:10.1128/JVI.01943-14.
- 18 A "Trojan horse" bispecific-antibody strategy for broad protection against ebolaviruses, Wec AZ, et al., *Science*, 354:350-354, 2016, doi:10.1126/science.aag3267.
- 19 Functional assessment of cell entry and receptor usage for SARS-CoV-2 and other lineage B betacoronaviruses, Letko M, Marzi A and Munster V, *Nat Microbiol*, 5:562-569, 2020, doi:10.1038/s41564-020-0688-y.
- 20 Structural basis for human coronavirus attachment to sialic acid receptors, Tortorici MA, et al., *Nat Struct Mol Biol*, 26:481-489, 2019, doi:10.1038/s41594-019-0233-y.
- 21 Potent MERS-CoV fusion inhibitory peptides identified from HR2 domain in spike protein of bat coronavirus HKU4, Xia S, et al., *Viruses*, 11, 2019, doi:10.3390/v11010056.
- 22 Towards a solution to MERS: protective human monoclonal antibodies targeting different domains and functions of the MERS-coronavirus spike glycoprotein, Widjaja I, et al., *Emerg Microbes Infect*, 8:516-530, 2019, doi:10.1080/22221751.2019.1597644.
- 23 Importance of neutralizing monoclonal antibodies targeting multiple antigenic sites on the Middle East Respiratory Syndrome coronavirus spike glycoprotein to avoid neutralization escape, Wang L, et al., *J Virol*, 92, 2018, doi:10.1128/JVI.02002-17.
- 24 Structural definition of a neutralization epitope on the N-terminal domain of MERS-CoV spike glycoprotein, Zhou H, et al., *Nat Commun*, 10:3068, 2019, doi:10.1038/s41467-019-10897-4.
- 25 Ultrapotent human neutralizing antibody repertoires against Middle East Respiratory Syndrome coronavirus from a recovered patient, Niu P, et al., *J Infect Dis*, 218:1249-1260, 2018, doi:10.1093/infdis/jiy311.
- 26 Evaluation of candidate vaccine approaches for MERS-CoV, Wang L, et al., *Nat Commun*, 6:7712, 2015, doi:10.1038/ncomms8712.
- 27 A humanized neutralizing antibody against MERS-CoV targeting the receptor-binding domain of the spike protein, Li Y, et al., *Cell Res*, 25:1237-1249, 2015, doi:10.1038/cr.2015.113.



ForteBio
47661 Fremont Boulevard
Fremont, CA 94538
888.OCTET-75 or 650.322.1360
ForteBio_info@sartorius.com

**ForteBio Analytics
(Shanghai) Co., Ltd.**
No. 88 Shang Ke Road
Zhangjiang Hi-tech Park
Shanghai, China 201210
fortebiosh@sartorius.com

Sartorius UK Ltd.
Longmead Business
CentreBlenheim Road, Epsom
Surrey KT19 9QQ
+44.1372.737100
Fortebio_europe@sartorius.com

Sartorius AG
Otto-Brenner-Straße 20
Göttingen 37079
+49.551.308.0
Fortebio_europe@sartorius.com