



Comparative data:

Sistemas de documentação de gel básicos:

Applications	InGenius	U:Genius3	G:BOX	NuGenius
Fluorescence UV transilluminator	Option	Option	Option	Option
Fluorescence UV Epi			Option	
Blue light	Option	Option	Option	Option
Visible stain converter	Option	Option	Option	Option
Colony plate	Option	Option	Option	Option
Spot blot	Option	Option	Option	
TLC			Option	

Sistemas de documentação de gel avançados:

Applications	GeneGnome XRQ	G:BOX Chemi XRQ	G:BOX Chemi XX6	G:BOX Chemi XX9	PXi/PXi Touch
Fluorescence UV transilluminator		Option	Option	Option	Option
Fluorescence UV Epi		Option	Option	Option	Option
Blue light		Option	Option	Option	Option
Visible stain converter		Option	Option	Option	Option
Colony plate		Option	Option	Option	Option
Spot blot		Option	Option	Option	Option
Fluorescence (colour)		Option	Option	Option	Option
Chemiluminescence	✓	✓	✓	✓	✓
Fluorescence IR			Option	Option	Option
Proteomics (2D)			Option	Option	Option



Paralab – Equi. Industriais e de Laboratório
www.paralab-bio.pt

email: info@paralab-bio.pt

Porto: Travessa do Calvário, 65, Giesta
 4420-392 - Valbom

Tel.: 224 664 320 Fax: 224 664 321