

Specifications

General

Sample volume:

50-120 ml

Particle size:

Up to 600µm – sample dependent

Stirrer speed:

350 - 3,500rpm

Dimensions (W, D, H):

Controller unit: 225mm x 170mm x 80mm Dispersion unit: 140mm x 175mm x 390mm

Weight:

Controller unit: 1kg Dispersion unit: 6.75kg

Power:

100-240V AC 50-60Hz 20W

Chemical compatibility:

Perlast® seals as standard. Full chemical compatibility requires the measurement cell to be fitted with an additional Perlast® seal.