

Specifications

Calibrant 1: Hi-Q Micro sensor calibrant

Standard type:

Duke 4010A.

Mean Size:

1.034 μ m.

Density:

1.050g/cc.

Target concentration:

9.0 x 10 particles/mL.

Volume:

20mL.

Calibrant 2: Hi-Q Nano sensor calibrant

Standard type:

Duke 3520A.

Mean Size:

0.519 μ m.

Density:

1.050g/cc.

Target concentration:

2.0 x 10 particles/mL.

Volume:

20mL.