

Processing time	45-60 min
Processing capability	1–12 samples per run
Extraction technology	magnetic particle separation technology
Protocol	programmed by scanning a barcode
Protocol input	barcode scanner
Sample volume	10–2000 μ l (depending on the kit)
Elution volume	50–400 μ l
Connection to PC	not required
Display	LCD (20 \times 4)
Certification	CE IVD
Nominal operating voltage	110–240 V, 50/60 Hz
Dimensions (W \times D \times H)	560 \times 590 \times 510 mm
Weight	55 kg