


CE

# C-Chip<sub>4ch</sub>

*Count your cells,  
for the most affordable price!*



# C-Chip<sub>4ch</sub>



Scan me



- ✓ New 4channel type disposable hemocytometer
- ✓ Created by the best-known microfluidic slide manufacturer
- ✓ High Quality, Low Price

## C-Chip 4channel is a Disposable Plastic Hemocytometer

- No need for coverslips
- Disposable, eliminates the need for washing
- Reduces exposure to infectious samples and risk of contamination
- Guaranteed reproducibility and reliability
- Four tests
- Made of sturdy quartz grade optical plastic

## Models

No.	Grid	Depth
DHC-N04	Neubauer Improved	100 $\mu$ m
DHC-B04	Burker Turk	100 $\mu$ m
DHC-F04	Fuchs Rosenthal	200 $\mu$ m




## Applications

Blood analysis (Hematology) : Blood cell counting  
Cell culture : Cell concentration measurement / Cell viability test  
Microbiology : Fungal spore counting

## Specifications

Dimension : 25 (W) x 75 (L) x Thickness 1.6 mm  
Chamber depth : 100  $\mu$ m, 200  $\mu$ m  
Chamber volume : 10  $\mu$ L (DHC-N04, B04)  
20  $\mu$ L (DHC-F04)



[website](http://www.nanoentek.com) www.nanoentek.com  
[e-mail](mailto:sales@nanoentek.com) sales@nanoentek.com

NESCT-N04-001E (V.1.3)

### NanoEntek, Inc.

Head Office  
12F, 5, Digital-ro 26-gil, Guro-gu, Seoul, 08389, Korea  
Tel +82-2-6220-7940 / Fax +82-2-6220-7999

### NanoEntek America, Inc.

220 Bear Hill Road, Suite 102, Waltham, MA 02451, USA  
Tel +1-781-472-2558 / Fax + 1-781-790-5649

NanoEntek Europe | mts med-tech supplies GmbH  
Lochhamerstr. 4a, 82152 Martinsried, Germany  
Tel +49-89-21-55-38-43 / Fax +49-89-99-95-46-60